
LUMBAR SUPPORT MASSAGER

QS300

USER MANUAL

MULTILINGUAL MANUALS

EN

Scan QR code for supported language manuals (Google recommended)
Scanne le code QR pour obtenir les manuels dans les langues prises en charge (recommandé par Google)
QR-Code scannen für Anleitungen in unterstützten Sprachen (Google empfohlen)
Scansion il codice QR per i manuali nelle lingue supportate (si consiglia Google)
Escane el código QR para manuales en idiomas admitidos (se recomienda Google)
Scan QR-koden for at få manualer på understøttede sprog (Google anbefales)
Skanna QR-koden för manualer på stödta språk (Google rekommenderas)
Naskenjuže QR kód pro příručky v podporovaných jazycích (Doporučujeme Google)
Σκάναρ το κωδικό QR για εγχειρίδια στις υποστηρίξιμες γλώσσες (συνιστάται η χρήση του Google)
Imbas kod QR untuk manual dalam bahasa yang disokong (dijayakan Google)

PRECAUTIONS

EN

To reduce the risk of interference and potential harm, do not use this product in conjunction with the following medical electronic devices:

- Implanted medical electronic devices, such as pacemakers.
- Life-sustaining medical electronic devices, such as artificial hearts and lungs.
- Medical electronic monitoring devices, such as electrocardiograms.

Prior to use, please ensure the main host buttons are functioning properly. If you are using the product for the first time or after the battery has been depleted from long-term storage, charge the product according to the instruction manual to ensure it works properly. Please dispose of this product in line with your city's environmental protection regulations to prevent environmental pollution.

- Do not use this product while driving, climbing, or engaging in other dangerous activities, as it may lead to accidents.
- Keep this product out of reach of infants, children, or individuals who are entirely or partially unconscious. Do not let them use it, as it may lead to accidents.
- This product is recommended for use by individuals aged 16+.
- Gradually increase the massage intensity according to personal comfort; otherwise, the discomfort from a higher intensity might cause muscle strain.
- This is a personal care massage product and cannot be used to substitute medicines or medical devices.
- If discomfort occurs during use, please stop using it immediately.
- If you have a waist injury, consult a healthcare professional before using this product.

In the following situations, please consult a physician before use; otherwise, accidents or discomfort may occur:

- Patients with infectious diseases.
- Patients allergic to stainless steel electrode pads.
- Patients with cardiac or cranial nerve disorders.
- Patients with sensory impairments or skin abnormalities.
- Patients who are unable to express their own will freely.
- Patients with perception disorders due to severe peripheral circulation problems caused by conditions such as diabetes.
- Patients with sensitive skin or allergic dermatitis.
- Individuals with circulatory disorders such as cerebrovascular diseases, inflamed or swollen waist skin, suppurative, or waist-related ailments should not use this product, as it may exacerbate their symptoms.
- For users who are pregnant, it is best to consult a doctor before use. Pregnant women should avoid using this product without the express permission of a doctor.

WARNING: Do not use this appliance near bathtubs, showers or water-filled reservoirs

OPERATION GUIDE

EN

Power button:
Long press for 1.5 seconds to power on (with voice prompt, indicator light on); long press again for 1.5 seconds to power off (with voice prompt, indicator light off). After powering on, short press to cycle through modes with voice prompts.

Hot compress button:
Short press to turn on/off hot compress and switch heat levels (default medium level). Long press for 1.5 seconds to enable/disable voice prompts (voice announces "Voice mode on"/"Voice mode off").

Working modes:
Relax Mode Revive Mode Hot Compress Mode (low/medium/high) De-stress Mode

BATTERY INDICATOR

EN

No.	Light Status	Description
1	Four lights steady on	Battery 75-100%
2	Three lights steady on	Battery 50-75%
3	Two lights steady on	Battery 25-50%
4	One light flashes	Battery 10-25%
5	One light flashes 3 times then shuts down	Battery 0-9%
6	Breathing light flashes steadily	Charging in progress
7	All lights steady on	Fully charged
8	One light flashes 5 times then shuts down	Power off for abnormal temperature

PRODUCT SPECIFICATIONS

EN

PRODUCT NAME: LUMBAR SUPPORT MASSAGER	
Model: QS300	Rated Input Voltage: 5V ⎓
Battery Capacity: 4500mAh/16.65Wh	Rated Input Power: 10W

PRECAUTIONS FOR PRODUCT DISPOSAL

EN

Please support the recycling of rechargeable battery (lithium-ion battery).

Danger

- DO NOT throw it into a fire, or heat it up to high temperature.
- DO NOT pierce it with nails, impact it by gravity or disassemble it for modification.
- DO NOT connect the positive and negative poles of the battery directly with metal.
- DO NOT store it together with necklaces, hairpins, etc.
- DO NOT charge, use or place the battery in areas with high temperatures, such as near a fireplace or under direct sunlight.
- The above situations may cause heating, fire and breakage.

Warning

Keep the removed rechargeable batteries out of reach of children to prevent accidental ingestion.

- If a battery is swallowed, consult a doctor immediately.
- In the event of a battery fluid leak, avoid touching the fluid and handle according to safety procedures.
- If battery electrolyte enters the eyes, rinse with clean water and seek medical advice without rubbing your eyes.
- If battery electrolyte contacts the skin, rinse the area with water and get medical help if needed.

Only disassemble scrapped products to prevent accidents such as fires or electric shocks.

For battery repair, contact the dealer where you purchased the product. This does not include battery replacement.

The lithium-ion battery in this product is recyclable and should not be placed, surrendered, it to a professional recycling agency when the product is discarded. Always turn off the device before removing the battery for disposal.

To safely dispose of the battery, follow these simplified steps:

- Use specialized tools to detach the device's casing.
- Remove the internal circuit board and motor.
- Extract the battery, sever its wires with pliers, and remove it completely.

CERTIFICATION STATEMENT

EN

Note: This product has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when used in a residential environment. This product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a specific installation. If this product does cause harmful interference to radio or television reception, which can be determined by turning the product off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the distance between the device and the receiver.
- Connect the device to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device in compliance with the essential requirements and other relevant provisions of Directive (EU)2023/988. CAUTION: THERE IS A RISK OF EXPLOSION IF THE BATTERY IS REPLACED WITH AN INCORRECT TYPE. DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS.

EU Declaration of Conformity: SKG Health Technologies Co., Ltd. hereby declare that this device complies with the essential requirements and other relevant provisions of Directive (EU)2023/988. A copy of the EU Declaration of Conformity is available online at <https://www.skg.com/pages/skg-declaration-of-conformity>.

Old electrical appliances must not be disposed of with the residual waste; they must be disposed of separately. Disposal at public collection points is free for private individuals. The owner of old devices is responsible for bringing them to these collecting points or similar facilities. By making this small personal effort, you contribute to the recycling of valuable raw materials and the proper treatment of toxic substances. available online at <https://www.skg.com/pages/skg-declaration-of-conformity>.

Read operator's manual


Do not use this device in a bathtub, shower or water-filled reservoir

UK REP: **HEALTHY CONSULTING LTD**
Office: 101-102 The Quadrant Street, London, United Kingdom, EC2A 4PU
Email: info@healthyconsulting.co.uk
Website: www.healthyconsulting.co.uk

EU REP: **ACCESS COMMERCE SL**
CALLE RIO TORRES N.º 1, PLANTA 1, SEDESPA, 03014 S. CARLOS DE LA RIBERA, SPAIN
Email: info@accesscommerce.es
Website: www.accesscommerce.es

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